
Passy-Muir Special Event Webinar

Swallowing Series

Passy-Muir® Valve Use with the Head and Neck Cancer Population



Linda Stachowiak MS/CCCSLP BRS-S MD Anderson Cancer Center Orlando



Head and Neck Cancer

• Refers to a collection of cancers arising from a variety of sites and grouped together under the following categories:

Oral cavity:

- $\bullet Lips$
- •Buccal mucosa
- •Anterior tongue
- •Floor of the mouth
- •Hard palate
- •Upper gingiva •Lower gingiva
- •Retromolar trigone
- Pharynx: Larynx, : Oropharynx •Supraglottic larynx
- Nasopharynx
- •Hypopharynx
- vocal cords and the mucosa of the
- anterior and posterior commissures)

•Glottic larynx (true

Subglottic

Head and Neck Cancer

Nasal cavity and the paranasal

- <u>sinuses:</u>
- Maxillary
- •Ethmoid
- Sphenoid
- •Frontal sinuses

Major salivary glands:

- •Paired parotids
- •Submandibular
- •Sublingual glands
- •Minor salivary glands

Head and Neck Cancer

- Cancer of the head and neck can significantly affect a person's ability to:
 - Communicate
 - Breathe
 - Swallow

_
_
_
_

Head and Neck Cancer Staging

- The tumor, node, metastasis (TNM) staging system of the American Joint Committee on Cancer (AJCC) is used to classify cancers of the head and neck
 - T refers to the primary tumor T1-T4, Tx, Tis
 - N stands for regional and local lymph nodes N1-N3
 - M refers to distant metastasis Mo vs. M1

Head and Neck Cancer Staging

- The staging of the tumor is one mechanism that drives the treatment. Treatment can include one or more of the following options:
 - Surgery
 - · Radiation therapy
 - Chemotherapy


Head and Neck Cancer Treatment: Factors influencing the treatment

- Location and its proximity to critical structures (functional deficits)
- Co-morbidities, candidacy for various tx options
- Previous treatment(s)
- Cost
- Patient preference
- Availability of support staff as needed
- Compliance
- Histology of the cancer
- Ability to provide the treatment in a particular facility

Challenges with Use of the Passy-Muir® Valve with the Head and Neck Cancer Population


- Timing of the need for tracheostomy is critical to Passy-Muir® Valve use with this population
- When is the tracheostomy necessary?
 - Prior to treatment
 - During the treatment
 - After the treatment

Tracheostomy Prior to Treatment • When a head and neck cancer patient requires a tracheostomy prior to treatment, it is often due to a glottic or supraglottic tumor which limits airflow, necessitating an extra breathing passage. • Passy-Muir®Valve is often contraindicated



Laryngectomy patients cannot use a Passy-Muir® Valve!

• Due to the anatomical changes associated with total laryngectomy, the Passy-Muir® Valve can not be used on this population.





Images Courtesy of InHealth Technologies

Tracheostomy *During* Treatment

- Some large reconstructive surgeries result in significant edema of the oral cavity requiring tracheostomy for an additional airway until edema resolves or dissipates
- Passy-Muir® Valve may be contraindicated initially, but frequently recommended after edema dissipates




Tracheostomy **During** Treatment


- Treatment with various forms of chemotherapy and/or radiation can significantly affect the airway in both positive and negative ways during the course of treatment.
- Induction chemotherapy may be initiated as the first line of treatment to shrink large glottic/supraglottic lesions. Dramatic improvements in the airway shortly after the onset of induction chemotherapy can often be seen with large tumors.







Concurrent chemoradiation may contribute to laryngeal edema at various points of treatment and the Passy-Muir®Valve may be contraindicated.










Treatment Plan

- Start with induction chemotherapy with TPF
- Follow with concurrent chemoradiation
- Can continue with oral intake of a regular diet at this time
- SLP to reassess swallow function after induction chemo therapy
- Patient given swallowing exercises prior to the onset of concurrent chemoradiation
- Patient to be seen for re-evaluation approximatley half way through concurrent chemoradiation




Pre-treatment exercises

- The goal of pre-treatment exercises is to preserve swallowing function and maximize functional outcomes at the completion of treatment.
- Specifically, the exercises target the base of tongue, the pharyngeal musculature, airway protection and laryngeal elevation.
- Patients are advised to practice each exercise in a series
 of 5 repetitions, three times a day as long they are able
 to. Realistically, certain early side effects may prohibit
 follow-through as instructed.

Post-induction chemotreatment videostroboscopy Results: L TVC fully mobile Slight improvement in R TVC Improved glottic closure No obvious mass effect Significantly improved glottic airway Consider capping trials, no trach removal (has upcoming concurrent chemoradiation)

SLP treatment during concurrent chemoradiation

- Continual monitoring of Passy-Muir® Valve use and swallowing functioning through treatment.
- The goal is to maintain nutritional status during treatment. The speech pathologist's goal would be to keep the nutrition oral versus feeding tube placement.
- This patient was able to maintain all of his nutrition PO throughout treatment but required diet modifications and postural strategies to remain safe and free from developing an aspiration pneumonia.



Conclusion

- The Passy-Muir® Valve can play an important role in treatment of the Head and Neck cancer patient by allowing verbal communication and assisting in swallow function
- The Head and Neck Cancer patient can fluctuate in communication and swallowing function at various points in the treatment process and subsequently requires constant monitoring and alterations in the treatment plan
- Ultimately this results in optimizing quality of life!

Questions and Conclusion

- Thank you for attending the webinar.
- Please complete your course evaluation for CEU credit.

For additional questions, email:

• Linda Stachowiak lstachowiak@cfl.rr.com

Please Participate! Quick & Important Research Survey

- Be part of <u>understanding</u> & working on <u>standardizing</u> dysphagia therapy in head & neck cancer patients.
- To do this, Dr. Susan Langmore and her team at Boston University Medical Center ask you to <u>please</u> take a 5-10 minute multiple choice survey.
- To take the survey, e-mail the survey coordinator william.sokoloff@bmc.org who will send you a link to the online survey. Your responses are truly invaluable!

Thank you!!!